



RADIATION SAFETY

# APRON CATALOGUE



# 01

---

Lead Garments are worn by healthcare professionals to protect them from the harmful effects of x-ray radiation during medical & other procedures requiring the use of x-ray or fluoroscopy. Lead aprons are usually worn in combination with other products to provide complete protection from x-ray.

---



# APRONS

## QUALITY. SERVICE. INTEGRITY.

Proguard aprons are proudly hand-stitched in the USA with one collective mission: perfection. Our light-weight designs and ergonomic styles are the result of decades of manufacturing experience. Seamless finishes, clean stitching, and elegant curves are what make our aprons best-in-industry.

Fulfilling our commitment to high quality, our Proguard Aprons are made in the USA from locally sourced materials whenever possible. We follow the strict guidelines of the ISO Quality Management System and is approved to sell in major markets throughout the world.



# APRON ORDERING GUIDE

Frontal Protection Style Example  
Flexback Apron

Full Protection Style Example  
Vest Skirt Apron

## 1 STEP 1: CHOOSE YOUR APRON STYLE

*We offer two types of aprons: Frontal and Full Protection.*

As their names suggest, a Frontal apron only protects the front of your body, while a Full Protection apron offers both front, lateral, and back protection.

For a **lighter weight apron**, choose a Frontal style. Our most popular frontal apron is the Flexback apron but you can view all of our styles beginning on page 24.

For **greater protection**, choose a Full Protection apron. We recommend the Proguard Balance Vest and Skirt, which offers a balance of protection and comfort.

Other important factors that should be considered in choosing your apron style: types of surgery, duration, body type and personal needs.







## FRONTAL APRONS

### THE GOOD

- Great flexibility and mobility
- Much lighter than full protection aprons
- Best for shorter procedures
- Forgiving on size

### THE NOT-SO-GOOD

- No back protection
- Reduced lateral protection



## VEST SKIRT APRONS

### THE GOOD

- Excellent protection - front, back and lateral protection
- Best for long procedures with high exposure
- Better than frontals for sitting procedures

### THE NOT-SO-GOOD

- One of the heaviest apron styles
- Not forgiving on size



## WRAPAROUND APRONS

### THE GOOD

- Excellent protection - front, back and lateral protection
- Best for long procedures with high exposure
- Easier to don and promotes good posture

### THE NOT-SO-GOOD

- Heavier than frontal style aprons
- Not forgiving on size



# 2

## STEP 2: CHOOSE YOUR CORE MATERIAL

Each core material varies in composition, weight and level of protection. Available Core Materials are listed below. For a more in-depth look, please refer to our core materials section



Pb

### STANDARD LEAD

- **Features:** Affordable & Highly Protective
- **Composition:** Lead
- **Certifications:** ASTM & IEC Certified



### PROLITE

- **Features:** Highly Flexible, 16% lighter than lead
- **Composition:** Cerium and/or Rare Earths
- **Certifications:** ASTM & IEC Certified

# 3

## STEP 3: CHOOSE YOUR CORE MATERIAL THICKNESS

Higher thicknesses result in more protection but also increase weight. **Our standard thickness options are 0.25mm, 0.35mm, and 0.50mm LE with .50mm\* and .35mm\* LE being the most popular options.** We offer additional options upon request such as .75mm, 1.00mm, etc.

Standard thicknesses in the front (frontal aprons, vest, skirts, one-piece wraparounds) are 0.50mm LE, and standard back protection thicknesses (vest, skirt, one-piece wraparounds) are 0.25mm LE. Standard thicknesses for other accessory styles may vary.

\*Please be advised that some state or country requirements may be specific about what thicknesses are permitted.



4

## STEP 4: CHOOSE YOUR FABRIC COLORS

You will need to choose colors for the exterior, interior and binding thread of your apron. Need a fabric that performs? Consider our anti-microbial, Satinguard fabrics for an easy-to-clean option. You can view our fabric color options on pages 59-64.

*If you have your own fabric that you'd like to use, we can accommodate! Just take a look at this Louis Vuitton apron; a special request by one of our customers.*



5

## STEP 5: CHOOSE CUSTOMIZATIONS

Depending on your apron style, you'll have a variety of customization options to choose from. Some examples: Embroidery, axillary wings, logos, extra pockets, etc.


**Please note:** Essentials aprons are only available in standard sizes and cannot be customized. If you need a made-to-measure Vest Skirt or Wraparound apron, please select the Balance configuration.



# CORE MATERIALS

## K-EDGE TECHNOLOGY

---



Knowing what's inside your apron is as important as how that apron claims to protect you.

Each Proguard core material is carefully formulated to more efficiently attenuate portions of the photon spectrum and limit “gaps” in protection due to over-reliance on individual x-ray protective elements. We utilize proper elemental balance to address k-edge concerns for key diagnostic imaging ranges. Our advanced elastomers exhibit the best balance of toughness, flexibility, durability and cracking resistance.


- **Antimony/barium provide more efficient attenuation of the photon spectrum below the K-edge window of lead ( 35 - 88 keV)**
- **Higher Z elements, tungsten, bismuth and lead are more efficient for stopping higher energy radiation (>9/88 keV), and also cover the K-edge fluorescence window of antimony ( < 35 keV).**





## STANDARD LEAD

- **Weight:** 3400 gsm (0.25mm)
- **Composition:** Lead Plasticizers
- **Certifications:** ASTM & IEC Certified




AVERAGE LEAD APRON WEIGHTS 0.50MM LE (LBS)

|                 | Frontal | Vest & Skirt | Wraparound |
|-----------------|---------|--------------|------------|
| <b>X-Small</b>  | 7.17    | 11.16        | 12.44      |
| <b>Small</b>    | 7.68    | 12.74        | 14.40      |
| <b>Medium</b>   | 8.63    | 14.75        | 15.95      |
| <b>Large</b>    | 10.35   | 17.96        | 17.52      |
| <b>X-Large</b>  | 11.74   | 19.56        | 19.08      |
| <b>2X-Large</b> | 12.28   | 21.46        | 22.00      |
| <b>3X-Large</b> | 14.21   | 23.84        | -          |



## PROLITE

- **Weight:** 2850 gsm (0.25mm) 16% lighter than lead
- **Composition:** Cerium and/or Rare earths, Plasticizers
- **Certifications:** ASTM & IEC Certified



AVERAGE PROLITE APRON WEIGHTS 0.50MM LE (LBS)

|                 | Frontal | Vest & Skirt | Wraparound |
|-----------------|---------|--------------|------------|
| <b>X-Small</b>  | 5.95    | 9.32         | 10.33      |
| <b>Small</b>    | 6.37    | 10.64        | 11.95      |
| <b>Medium</b>   | 7.16    | 12.32        | 13.23      |
| <b>Large</b>    | 8.59    | 15.00        | 14.54      |
| <b>X-Large</b>  | 9.74    | 16.33        | 15.84      |
| <b>2X-Large</b> | 10.20   | 17.92        | 18.26      |
| <b>3X-Large</b> | 11.79   | 19.91        | -          |

# TAKING CARE OF YOUR APRON

## DO'S

Clean most dirt and stains (blood, bodily fluids, barium contrast media) using designated apron cleaning solutions or mild household soaps in diluted form. To apply cleaning fluid, please use a cloth or wipe and use room temperature water.

## DONT'S


While it is possible, we DO NOT recommend you machine wash or dry your apron. Do not use hot or warm water. Do not iron or autoclave.

## STORAGE

When not in use, garments should be placed on a hanger or stored using Proguard apron storage equipment. Avoid contact with hot surfaces or exposure to direct sunlight. Also avoid contact with sharp objects, as these may damage the garment's protective core material. For transportation, garments should be carefully rolled up, not folded, and placed into a box or Proguard duffle bag. All of our apron storage options can be found at the end of the apron section.

## LEAD DISPOSAL & RECYCLING

Aprons containing lead can be extremely hazardous to the environment and must be disposed of properly. We provide an in-house replacement service for this very occasion. Our disposal service is free with the purchase of any new apron! Simply call us or your local representative.





# 2 YEAR WARRANTY

Proguard high quality garments (“Garments”, aka PPE or Aprons) include; aprons, thyroid collars, gloves, dental aprons, reproductive aprons, and lead accessories. The Garments are warranted to be free of defects in materials and workmanship to the original purchaser for 2 years.

Upon evaluation, we will either repair or replace the defective item or part without charge to the purchaser. This warranty is void when the product has been tampered with, when repairs or attempted repairs have been made by unauthorized persons, or when the item has been subject to misuse, abuse or damage in transit. No charges will be accepted or reimbursed for unauthorized parts, repairs or services.

Please contact your local representative for return instructions. All returns, repairs, or exchanges, regardless of reason, must be accompanied by an authorized RMA. Customer should make available the Purchase Order # or relevant sales information.



# COMPLIANCE

**EN 61331-1:2014: (EU)** Protective devices against Diagnostic Medical X-ray radiation. (Determination of attenuation properties of material)

**EN 61331-3:2014: (EU)** Protective devices against diagnostic medical X-radiation Part 3: Protective clothing, Eyewear and Protective Patient shields.

**ASTM 2547-18: (USA)** Standard test method for determining the attenuation properties in a primary x-ray beam of materials used against radiation.



# MANUFACTURER

## FDA LISTING

OWNER/OPERATOR #: 10025412  
FDA LISTING NUMBER: D048258  
PRODUCT CODE: IWO

DEVICE CLASS: #1  
REGULATION #: 892.6500  
OSHA REGULATION: 1910.1096

## QUALITY SYSTEM

ISO 9001  
EU 2016/425, MODULE D, B

# #1 IN CUSTOM APRONS

Proguard specializes in tailor-made aprons. You can choose from over 100 colors for your apron's exterior color, interior color, and binding. When you're finished, you can add embroidery, logos and other add-ons for a truly customized experience.

Hand-crafted and made with patented polymer technology, our Proguard Aprons are more light-weight, flexible and durable than the competition.



PROGUARD







# PROGUARD



## A COMMITMENT TO QUALITY

With full time quality control staff and ISO certification, you are getting the best in American made protective apparel.



## ULTRA-FLEXIBLE CORE MATERIALS

Proguard's Core materials feature more flexibility and durability which provides superior, long lasting comfort.



## ECO-FRIENDLY LEAD FREE

Proguard lead free core materials are recyclable and non-toxic.



## LUXRAY PATTERNS

LuXRay offers patterned accents on the apron's sides and shoulders combined with an optional matching design on your thyroid collar.



## IN-HOUSE EMBROIDERY & LOGOS

Don't let your apron "disappear." Our embroidery and logos are done in-house with modern monogram machines and the latest in digitization technology.



## NON-SLIP FABRIC ON SKIRT INTERIOR

Included for free, our premium non-slip fabric improves the fit of your skirt by reducing sliding and shifting.

# CUSTOMIZATION OPTIONS



## BACK STRAP

The Velcro backstrap keeps your apron shoulders secured and prevents slippage.



## ELASTIC BACK

Provides enhanced adjustability and lumbar support by shifting apron weight from shoulders to the hips.



## BACK MESH

Allows for maximum airflow in the back panels of your frontal apron and helps keep you cool.



## NON-SLIP FABRIC

A high friction, anti-slip fabric that is added to the back interior of your skirt or other key areas depending on style selection which reduces slippage and improves fit.





# CUSTOMIZATIONS: COMFORT OPTIONS

When it comes to designing your perfect custom apron, Proguard is #1. Our great selection of apron add-ons can increase your protection and improve comfort. We offer 3 types of apron customizations: Comfort, Protection and Design add-ons. Each customization is clearly defined to ensure you make the right choices for you.



## APRON SUSPENDERS



Need a little extra help holding up your skirt? Try our apron suspenders to keep that apron skirt right where you want it.



## SURGICAL DROP SHOULDERS

Specialty velcro tabs can be added to different apron styles allowing for easy removal of your apron to maintain sterile field.





# CUSTOMIZATIONS: DESIGN OPTIONS



## LOGO AND EMBROIDERY

Personalize your apron by placing your name or logo anywhere! Common locations include the pocket, just above the pocket, and on your thyroid.

## VELCRO PATCH

A customized, embroidered velcro patch can be added to a specific area on your apron. Typical locations are on the pocket, above pocket, or to the left of your pocket.



## LUXRAY PATTERNS

Have your apron tailored to perfection with LuxRay's patterned accents on the apron sides, pocket and thyroid collar. If you have a special fabric you'd like to use, we can incorporate that as well!

## EXTRA POCKETS

Need an extra pocket or two? Whether you need one on your skirt, inside your vest or even on your knees (yes, we did it), we can accommodate most requests. Shown above is a spacious L-shaped pocket which can be added to your skirt.



# CUSTOMIZATIONS: PROTECTION OPTIONS



## AXILLARY WINGS


---

Axillary wings are a protective flap sewn-in or “snapped-on” to your apron and offer vital protection to vulnerable breast tissue and red bone marrow areas. Protection levels will match your apron specifications.

## SIDE EXTENSIONS

---

Side Extensions can provide an additional +2 inches of protection on both sides and provide important lateral protection. Protection levels will match your apron specifications.



## APRON SLEEVE

---

Apron sleeves wrap around your entire arm and offer comprehensive protection of vulnerable breast tissue and red bone marrow areas. Sleeves come standard in 0.25mm LE and can be sewn-in permanently or offered in our removable “Snap-on” option.

**Note:** To improve mobility, our apron sleeves are not sewn/attached at the armpit. For an upcharge, we can add a small “flap” which is tucked inside the apron armpit to ensure there are no gaps in protection.





APRON INSPIRATION

# CUSTOM PROJECTS



## LOUIS VUITTON FLEXBACK APRON

One of our favorite custom apron projects was this LuxRay Flexback featuring Louis Vuitton accents. Custom from top to bottom, this apron boasted an attached thyroid collar and velcro-adjustable apron sleeve also made in the Louis Vuitton fabric.

## SWISS SOCCER CLUB SKIRT

This Swiss soccer club apron required great attention to detail to complete. The customer provided a special team Flag to be utilized as the actual apron fabric.

If you have a special fabric you want to use, send it our way! Our team loves a challenge.







## THE DAWN OF JUSTICE APRONS

After custom flags and Louis Vuitton it was necessary that we honor the age of super heroes with our Special Edition DC Comics series. We hope Marvel and Avatar themes are next!





# PROGUARD



## FRONT PROTECTION

# APRON FEATURES



### Backstrap

Velcro backstrap securely fastens shoulder area.

### Shoulder Padding

Premium foam padding in the shoulders provides extra comfort.



### Back Relief Belt

Wide, elastic belt provides lumbar support and maximum comfort.

### Apron Pocket

### Belt Buckle Closure

## ADDITIONAL CUSTOMIZATION OPTIONS

|                |               |           |
|----------------|---------------|-----------|
| Thyroid Collar | Drop Shoulder | Wings     |
| Extra Pocket   | Embroidery    | Back Mesh |
| Logo           | Sleeve        |           |

## Core Material Weight (0.25mm) Thickness Options

| Core Material | Weight (0.25mm) | Thickness Options  |
|---------------|-----------------|--------------------|
| Pb Lead       | 3400 gsm        | 0.25, 0.50mm       |
| Prolite       | 2850 gsm        | 0.25, 0.35, 0.50mm |
| Premier       | 2600 gsm        | 0.25, 0.35, 0.50mm |
| Nano          | 2300 gsm        | 0.25, 0.35, 0.50mm |



FOR SIZING AND COLOUR OPTIONS SEE PAGES 55-60.



\*CE applies to Lead and Prolite in 0.25, 0.35, and 0.50mm. Prolite is limited to 110kVp diagnostic range.

# BACK RELIEF APRON

MODEL: [ AP-BR ]

Save yourself from back pain with the Proguard Back Relief apron. This apron features a unique padded belt system constructed with a bonded elasticised panel that will not lose its shape. When engaged, the special back relief belt provides lumbar support while transferring the weight of the apron from your shoulders to your hips for increased comfort and back support throughout the life of the apron. A velcro backstrap is also included to secure the apron around your shoulder area and prevent sliding.





# BUCKLE APRON

MODEL: [ AP-BA ]

The Buckle Apron comes standard with an easy-release buckle allowing for quick removal of the apron. The belt is adjustable to accommodate different sizes. Criss-cross back straps provide a secure, comfortable fit.



## FRONT PROTECTION

### APRON FEATURES



**Wide Body Back Panels**  
The comfortable panels secure your apron and provides lumbar support.

**Shoulder Padding**  
Premium foam padding in the shoulders provides extra comfort.

**Apron Pocket**

**Durable Buckle Closure**

### ADDITIONAL CUSTOMIZATION OPTIONS

|                |               |              |
|----------------|---------------|--------------|
| Thyroid Collar | Drop Shoulder | Elastic Back |
| Extra Pocket   | Embroidery    | Backstrap    |
| Logo           | Sleeve        |              |

| Core Material | Weight (0.25mm) | Thickness Options  |
|---------------|-----------------|--------------------|
| Lead          | 3400 gsm        | 0.25, 0.50mm       |
| Prolite       | 2850 gsm        | 0.25, 0.35, 0.50mm |



FOR SIZING AND COLOUR OPTIONS SEE PAGES 55-60.



\*CE applies to Lead and Prolite in 0.25, 0.35, and 0.50mm. Prolite is limited to 110kVp diagnostic range.

FRONT PROTECTION

# APRON FEATURES

**+** **Backstrap**  
Velcro backstrap securely fastens shoulder area.

**+** **Elastic Back**  
Provides enhanced adjustability and lumbar support. Shifts weight from shoulders & hips.

**Shoulder Padding**  
Premium foam padding in the shoulders provides extra comfort.

**Apron Pocket**

**Velcro Closure**

## ADDITIONAL CUSTOMIZATION OPTIONS

- |                |               |           |
|----------------|---------------|-----------|
| Thyroid Collar | Drop Shoulder | Wings     |
| Extra Pocket   | Embroidery    | Back Mesh |
| Logo           | Sleeve        |           |

| Core Material | Weight (0.25mm) | Thickness Options |
|---------------|-----------------|-------------------|
|---------------|-----------------|-------------------|

| Core Material | Weight (0.25mm) | Thickness Options  |
|---------------|-----------------|--------------------|
| Lead          | 3400 gsm        | 0.25, 0.50mm       |
| Prolite       | 2850 gsm        | 0.25, 0.35, 0.50mm |



FOR SIZING AND COLOUR OPTIONS SEE PAGES 55-60.



\*CE applies to Lead and Prolite in 0.25, 0.35, and 0.50mm. Prolite is limited to 110kVp diagnostic range.



# FLEXBACK APRON

MODEL: [ AP-FB ]

Representing the gold standard of quality, the Flexback apron design provides enhanced back relief, comfort, and excellent weight distribution. Criss-cross back straps are specially designed with elastic to “flex” the apron’s weight burden off your shoulders to your hips, thereby providing support to your lower back. A velcro backstrap (shown below) is included and secures the apron around your shoulder area to prevent sliding.



# SURGICAL DROP APRON

MODEL: [ AP-SD ]

Maintain your sterile field with the quick release Surgical Drop Apron. Easy-release velcro closures allow the surgical drop to fall even when underneath your scrubs. With padded shoulders and a criss-cross back strap design the surgical drop apron provides maximum comfort for all users. Hanging straps are included to allow for easy storage on any Proguard apron rack.



## FRONT PROTECTION

### APRON FEATURES

- +** **Hanging Straps**  
Hanging straps on the shoulders allow for easy storage on any Proguard apron rack.
- +** **Shoulder Padding**  
Premium foam padding in the shoulders provides extra comfort.
- +** **Surgical Drop**  
Quick release velcro closure tabs help maintain sterile field.
- +** **Apron Pocket**
- +** **Velcro Closure**

### ADDITIONAL CUSTOMIZATION OPTIONS

- Thyroid Collar
- Extra Pocket
- Embroidery
- Logo

### Core Material      Weight (0.25mm)      Thickness Options

|                  |          |                    |
|------------------|----------|--------------------|
| <b>Pb</b> Lead   | 3400 gsm | 0.25, 0.50mm       |
| <b>P</b> Prolite | 2850 gsm | 0.25, 0.35, 0.50mm |



FOR SIZING AND COLOR OPTIONS SEE PAGES 55-60.



\*CE applies to Lead and Prolite in 0.25, 0.35, and 0.50mm. Prolite is limited to 110kVp diagnostic range.



FRONT PROTECTION



# APRON FEATURES

- Shoulder Padding  
Premium foam padding in the shoulders provides extra comfort.
- Apron Pocket
- Tie Closure

## ADDITIONAL CUSTOMIZATION OPTIONS

- |                |            |           |
|----------------|------------|-----------|
| Thyroid Collar | Embroidery | Wings     |
| Logo           | Sleeve     | Backstrap |

| Core Material | Weight (0.25mm) | Thickness Options |
|---------------|-----------------|-------------------|
|---------------|-----------------|-------------------|

| Core Material   | Weight (0.25mm) | Thickness Options  |
|---|-----------------|--------------------|
|  Lead    | 3400 gsm        | 0.25, 0.50mm       |
|  Prolite | 2850 gsm        | 0.25, 0.35, 0.50mm |



FOR SIZING AND COLOR OPTIONS SEE PAGES 55-60.



\*CE applies to Lead and Prolite in 0.25, 0.35, and 0.50mm. Prolite is limited to 110kVp diagnostic range.

# TIE APRON

MODEL: [ AP-TA ]

The Tie Apron is an easy-to-use, affordable option for x-ray protection. Utilizing a tie closure, the Tie Apron's criss crossing straps ensure a secure and comfortable fit. See Pediatric options at the end of this section.



# UROLOGY APRON

MODEL: [ AP-UA ]

If you sit while performing procedures requiring x-ray, you may be unnecessarily exposing yourself to radiation. The Urology Apron features a three panel system at the bottom of the apron which enables the user to sit comfortably while maintaining maximum coverage from x-ray.



## FRONT PROTECTION

## APRON FEATURES

- +** **Multi-Panel Design**  
Allows for easy maneuvering and a comfortable fit.
- +** **Backstrap**  
Velcro backstrap securely fastens shoulder area.
- +** **Shoulder Padding**  
Premium foam padding in the shoulders provides extra comfort.
- +** **Apron Pocket**
- +** **Velcro Closure**

## ADDITIONAL CUSTOMIZATION OPTIONS

- Thyroid Collar
- Drop Shoulder
- Wings
- Extra Pocket
- Embroidery
- Back Mesh
- Logo
- Sleeve

| Core Material | Weight (0.25mm) | Thickness Options  |
|---------------|-----------------|--------------------|
| Lead          | 3400 gsm        | 0.25, 0.50mm       |
| Prolite       | 2850 gsm        | 0.25, 0.35, 0.50mm |



FOR SIZING AND COLOR OPTIONS SEE PAGES 55-60.



\*CE applies to Lead and Prolite in 0.25, 0.35, and 0.50mm. Prolite is limited to 110kVp diagnostic range.

FRONT PROTECTION

# APRON FEATURES



**Wide Body Back Panels**

The comfortable panels secure your apron and provides lumbar support.

**Shoulder Padding**

Premium foam padding in the shoulders provides extra comfort.

**Apron Pocket**

**Velcro Closure**

## ADDITIONAL CUSTOMIZATION OPTIONS

- |                |               |           |
|----------------|---------------|-----------|
| Thyroid Collar | Drop Shoulder | Wings     |
| Extra Pocket   | Embroidery    | Backstrap |
| Logo           | Sleeve        | Back Mesh |

| Core Material | Weight (0.25mm) | Thickness Options |
|---------------|-----------------|-------------------|
|---------------|-----------------|-------------------|

| Core Material | Weight (0.25mm) | Thickness Options  |
|---------------|-----------------|--------------------|
| Lead          | 3400 gsm        | 0.25, 0.50mm       |
| Prolite       | 2850 gsm        | 0.25, 0.35, 0.50mm |



FOR SIZING AND COLOR OPTIONS SEE PAGES 55-60.



\*CE applies to Lead and Prolite in 0.25, 0.35, and 0.50mm. Prolite is limited to 110kVp diagnostic range.

# VELCRO ADJUSTABLE APRON

MODEL: [ AP-VA ]

Do you need comfort and adjustability at an affordable price? The Velcro Adjustable apron has crisscrossing back straps that provide back support while securing the apron around your waist. Wide back straps prevent the apron from slipping around your shoulder and back area. Pediatric sizing available, see the end of the apron section.





# MATERNITY APRON

MODEL: [ AP-MA ]

The maternity apron features easy ties and a quick release buckle to secure the back and front panels; allowing complete adjustability during pregnancy. The maternity apron also includes an additional .50mm equivalent half apron for a total of 1.00mm Pb protection over the abdominal area and lower extremities.




FULL PROTECTION

## APRON FEATURES

- + Complete Coverage**  
The front & back panels ensure you are fully protected from beginning to end of your pregnancy.
- + Half Apron Included**  
A Half apron is included for additional protection as you progress through your pregnancy.
- Shoulder Padding**  
Premium foam padding in the shoulders provides extra comfort.
- Apron Pocket**
- Tie & Buckle Closures**

### ADDITIONAL CUSTOMIZATION OPTIONS

- Thyroid Collar
- Wings
- Extra Pocket
- Embroidery
- Logo
- Sleeve

| Core Material   | Weight (0.25mm) | Thickness Options  |
|---|-----------------|--------------------|
|  Prolite | 2850 gsm        | 0.25, 0.35, 0.50mm |






FOR SIZING AND COLOR OPTIONS SEE PAGES 55-60.



\*CE applies to Lead and Prolite in 0.25, 0.35, and 0.50mm. Prolite is limited to 110kVp diagnostic range.

FULL PROTECTION

## APRON FEATURES


-  **Elastic Back**  
Provides enhanced adjustability and lumbar support. Shifts weight from shoulders & hips.
-  **Single Front Panel**  
Eliminates excess overlapping panels achieving a 15% weight reduction (compared to other wraparound and VS models).
-  **Tri-tab Closure**  
Ensures a proper fit.

- Shoulder Padding**  
Premium foam padding in the shoulders provides extra comfort.
- Apron Pocket**
- Velcro Closure**

### ADDITIONAL CUSTOMIZATION OPTIONS

- |                |       |              |
|----------------|-------|--------------|
| Thyroid Collar | Wings | Extra Pocket |
| Embroidery     | Logo  | Sleeve       |

| Core Material | Weight (0.25mm) | Thickness Options |
|---------------|-----------------|-------------------|
|---------------|-----------------|-------------------|

|   |          |                    |
|---|----------|--------------------|
|  Prolite | 2850 gsm | 0.25, 0.35, 0.50mm |
|---|----------|--------------------|

# WRAPAROUND FLEXBACK MATERNITY APRON

MODEL: [ AP-WAFBMA ]

The Wraparound Flexback Maternity apron combines the best qualities of our flexback apron with the full protection and back support of our wraparound back relief apron. The elastic stretch-back shifts weight from your shoulders to your hips while the tri-tab velcro closure system provides convenient, easy-on easy-off wearability. In addition, the 6" wide frontal velcro allows the wearer to adjust the crossties up or down to account for mom's growing belly. The Wrap Around Flexback eliminates excess overlapping lead and is up to **15% lighter than traditional wrap styles.**



FOR SIZING AND COLOR OPTIONS SEE PAGES 55-60.



\*CE applies to Lead and Prolite in 0.25, 0.35, and 0.50mm. Prolite is limited to 110kVp diagnostic range.

# ESSENTIAL™ WRAPAROUND BACK RELIEF APRON




MODEL: [ AP-WAQ ]

Offering full protection, this back relieving design transfers the weight from your shoulders to your hips with the aid of a 6" lumbar support belt. The Wrap Around Back Relief apron features velcro and buckle closures to ensure a proper fit. The Essential apron configuration **features extended front panels to achieve a 3/4 overlap (where protection is doubled) which allows for a lighter weight apron without compromising protection. The Essential Wraparound is available in standard sizing only.** For a made-to-measure Wrap-around apron, see the Balance Wraparound on the next page.



**Example of a 0.50mm LE  
Essential Wraparound  
Apron Lead Configuration**



FULL PROTECTION

## APRON FEATURES

- +** **Adjustable Buckle Set**  
Adjustable buckle closures ensure a secure fit.
- +** **Back Relief Belt**  
Elastic lumbar support belt with buckle closure system.
- +** **Extended Front Panels**  
Front panels are extended to 3/4 overlap coverage where LE protection is doubled
- +** **Shoulder Padding**  
Premium foam padding in the shoulders provides extra comfort.
- +** **Apron Pocket**
- +** **Velcro and Buckle Closures**

### ADDITIONAL CUSTOMIZATION OPTIONS

- Thyroid Collar
- Wings
- Extra Pocket
- Embroidery
- Logo
- Sleeve

| Core Material | Weight (0.25mm) | Thickness Options  |
|---------------|-----------------|--------------------|
| Lead          | 3400 gsm        | 0.25, 0.50mm       |
| Prolite       | 2850 gsm        | 0.25, 0.35, 0.50mm |



FOR SIZING AND COLOR OPTIONS SEE PAGES 55-60.



\*CE applies to Lead and Prolite in 0.25, 0.35, and 0.50mm. Prolite is limited to 110kVp diagnostic range.



FULL PROTECTION

# APRON FEATURES

- +** **Adjustable Buckle Set**  
Adjustable buckle closures ensure a secure fit.
- +** **Back Relief Belt**  
Elastic lumbar support belt with buckle closure system.
- +** **Seam-to-seam Protection**  
Seam-to-seam, standardized protection in Front (no change in protection at overlap)
- +** **Shoulder Padding**  
Premium foam padding in the shoulders provides extra comfort.
- +** **Apron Pocket**
- +** **Velcro and Buckle Closures**

## ADDITIONAL CUSTOMIZATION OPTIONS

- Thyroid Collar
- Wings
- Extra Pocket
- Embroidery
- Logo
- Sleeve

| Core Material | Weight (0.25mm) | Thickness Options  |
|---------------|-----------------|--------------------|
| Lead          | 3400 gsm        | 0.25, 0.50mm       |
| Prolite       | 2850 gsm        | 0.25, 0.35, 0.50mm |



FOR SIZING AND COLOR OPTIONS SEE PAGES 55-60.



\*CE applies to Lead and Prolite in 0.25, 0.35, and 0.50mm. Prolite is limited to 110kVp diagnostic range.

# BALANCE™ WRAPAROUND BACK RELIEF APRON

MODEL: [ AP-WABR ]

Offering full protection, this back relieving design transfers the weight from your shoulders to your hips with the aid of a 6" lumbar support belt. The Wrap Around Back Relief apron features velcro and buckle closures to ensure a proper fit. The Balance configuration **features seam-to-seam, single LE protection with standard overlapping panels** (see below for lead configuration).

Example of a 0.50mm LE Balance Wraparound Apron Lead Configuration



# WRAPAROUND FLEXBACK APRON




MODEL: [ AP-WAFB ]

Want the ultimate in comfort and protection? The Wrap Around Flexback apron combines the best qualities of our flexback apron with the full protection and back support of our wraparound back relief apron. The elastic stretch-back shifts weight from your shoulders to your hips while the tri-tab velcro closure system provides convenient, easy-on easy-off wearability. The Wrap Around Flexback eliminates excess overlapping lead and is up to **15% lighter than traditional wrap styles.**





FULL PROTECTION

## APRON FEATURES

-  **Elastic Back**  
Provides enhanced adjustability and lumbar support. Shifts weight from shoulders & hips.
-  **Single Front Panel**  
Eliminates excess overlapping panels achieving a 15% weight reduction (compared to other wraparound and VS models).
-  **Tri-tab Closure**  
Ensures a proper fit.
- Shoulder Padding**  
Premium foam padding in the shoulders provides extra comfort.
- Apron Pocket**
- Velcro Closure**

### ADDITIONAL CUSTOMIZATION OPTIONS

- Thyroid Collar
- Wings
- Extra Pocket
- Embroidery
- Logo
- Sleeve

| Core Material   | Weight (0.25mm) | Thickness Options  |
|---|-----------------|--------------------|
|  Lead    | 3400 gsm        | 0.25, 0.50mm       |
|  Prolite | 2850 gsm        | 0.25, 0.35, 0.50mm |



FOR SIZING AND COLOR OPTIONS SEE PAGES 55-60.



\*CE applies to Lead and Prolite in 0.25, 0.35, and 0.50mm. Prolite is limited to 110kVp diagnostic range.

FULL PROTECTION

## APRON FEATURES

- + **Elastic Back**  
Provides enhanced adjustability and lumbar support. Shifts weight from shoulders & hips.
- + **Single Front Panel**  
Eliminates excess overlapping panels achieving a 10% weight reduction (compared to VS model).
- + **Tri-tab Closure Design**  
Ensures a proper fit.
- + **Non-slip Backing (skirt)**  
Reduces slippage and improves fit.
- Shoulder Padding**  
Premium foam padding in the shoulders provides extra comfort.
- Apron Pocket**
- Multi-tab Velcro Closure**
- Hanging Tabs on Skirt**

### ADDITIONAL CUSTOMIZATION OPTIONS

- |                |              |                    |
|----------------|--------------|--------------------|
| Thyroid Collar | Extra Pocket | Embroidery         |
| Logo           | Sleeve       | Suspenders (skirt) |

| Core Material | Weight (0.25mm) | Thickness Options |
|---------------|-----------------|-------------------|
|---------------|-----------------|-------------------|

|    |         |          |                    |
|----|---------|----------|--------------------|
| Pb | Lead    | 3400 gsm | 0.25, 0.50mm       |
|    | Prolite | 2850 gsm | 0.25, 0.35, 0.50mm |



FOR SIZING AND COLOR OPTIONS SEE PAGES 55-60.



\*CE applies to Lead and Prolite in 0.25, 0.35, and 0.50mm. Prolite is limited to 110kVp diagnostic range.



# TRI-TAB VEST AND SKIRT APRON

MODEL: [AP-TVEST]<sub>(vest only)</sub> [AP-SKIRT]<sub>(skirt only)</sub> [AP-TVS]<sub>(vest and skirt combo)</sub>

With an easy-on, easy-off multi tab (velcro) closure design, the Tri-tab is a popular style for vest and skirt combinations. Velcro closures are located on the left shoulder and around the hips. An elastic back allows for a flexible fit while providing premium comfort and back support. The Tri-tab Vest Skirt reduces excess overlapping lead and is **up to 10% lighter than traditional vest models.**






# ESSENTIAL™ VEST AND SKIRT APRON

MODEL: [AP-VQ]<sub>(vest only)</sub> [AP-SKIRT]<sub>(skirt only)</sub> [AP-VSQ]<sub>(vest and skirt combo)</sub>

Be prepared for prolonged procedures with the Vest Skirt apron. The two piece design shifts 70% of the weight directly to the hips, reducing back fatigue. The Essential apron configuration **features extended front panels to achieve a 3/4 overlap (where protection is doubled) which allows for a lighter weight apron without compromising protection.** The skirt comes standard with a premium non-slip fabric in the back. **Essential aprons are only available in standard sizes.** For a made-to-measure Vest Skirt, see the Balance Vest Skirt Apron.



**Example of a 0.35mm LE Essential Vest Skirt Apron Lead Configuration**



## FULL PROTECTION

### APRON FEATURES

- +** Adjustable Buckle Set  
Adjustable buckle closures ensure a secure fit.
- +** Non-slip Backing (skirt interior)  
Reduces slippage and improves fit.
- +** Extended Front Panels  
Front panels are extended to 3/4 overlap coverage where LE protection is doubled
- Shoulder Padding  
Premium foam padding in the shoulders provides extra comfort.
- Apron Pocket
- Velcro Closure
- Hanging Tabs on Skirt

### ADDITIONAL CUSTOMIZATION OPTIONS

- Thyroid Collar
- Wings
- Extra Pocket
- Suspenders (skirt)
- Embroidery
- Logo
- Sleeve

| Core Material | Weight (0.25mm) | Thickness Options  |
|---------------|-----------------|--------------------|
| Lead          | 3400 gsm        | 0.25, 0.50mm       |
| Prolite       | 2850 gsm        | 0.25, 0.35, 0.50mm |



FOR SIZING AND COLOR OPTIONS SEE PAGES 55-60.



\*CE applies to Lead and Prolite in 0.25, 0.35, and 0.50mm. Prolite is limited to 110kVp diagnostic range.

FULL PROTECTION



# BALANCE™ VEST AND SKIRT APRON

## APRON FEATURES

- +** **Adjustable Buckle Set**  
Adjustable buckle closures ensure a secure fit.
- +** **Non-slip Backing (skirt interior)**  
Reduces slippage and improves fit.
- +** **Seam-to-seam Protection**  
Seam-to-seam, standardized protection in Front (no change in protection at overlap)
- Shoulder Padding**  
Premium foam padding in the shoulders provides extra comfort.
- Apron Pocket**
- Velcro Closure**
- Hanging Tabs on Skirt**

### ADDITIONAL CUSTOMIZATION OPTIONS

- Thyroid Collar
- Wings
- Extra Pocket
- Suspenders (skirt)
- Embroidery
- Logo
- Sleeve

### Core Material      Weight (0.25mm)      Thickness Options

|  |         |          |                    |
|--|---------|----------|--------------------|
|  | Lead    | 3400 gsm | 0.25, 0.50mm       |
|  | Prolite | 2850 gsm | 0.25, 0.35, 0.50mm |



FOR SIZING AND COLOR OPTIONS SEE PAGES 55-60.




\*CE applies to Lead and Prolite in 0.25, 0.35, and 0.50mm. Prolite is limited to 110kVp diagnostic range.

MODEL: [AP-V]<sub>(vest only)</sub> [AP-SKIRT]<sub>(skirt only)</sub> [AP-VS]<sub>(vest and skirt combo)</sub>

Be prepared for prolonged procedures with the Vest Skirt apron. The two piece design shifts 70% of the weight directly to the hips, reducing back fatigue. The Balance configuration **features seam-to-seam, single LE protection with standard overlapping panels** (see below for lead configuration). The skirt comes standard with a premium non-slip fabric in the back.

### Example of a 0.35mm LE Balance Vest Skirt Apron Lead Configuration



# ESSENTIAL™ FLEXBACK VEST AND SKIRT APRON




MODEL: [AP-VFBQ]<sub>(vest only)</sub> [AP-SKIRT]<sub>(skirt only)</sub> [ AP-VSFBQ ]<sub>(vest and skirt combo)</sub>

The Essential Vest Skirt Flexback shifts 70% of the weight directly to the hips with added support of an adjustable, elastic back relief belt. The Essential apron configuration **features extended front panels to achieve a 3/4 overlap (where protection is doubled) which allows for a lighter weight apron without compromising protection.** The skirt comes standard with an anti-slip fabric in the back (interior) to keep your skirt in place. **Essential aprons are only available in standard sizes.** For a made-to-measure Vest Skirt Flexback, see the Balance Vest Skirt Flexback Apron.



**Example of a 0.50mm LE Essential Vest Skirt Apron Lead Configuration**



## FULL PROTECTION

### APRON FEATURES

- + **Adjustable Buckle Set**  
Adjustable buckle closures ensure a secure fit.
- + **Flexback Belt**  
Wraps internally around your back and securely fastens in the front.
- + **Non-slip Backing (skirt interior)**  
Reduces slippage and improves fit.
- + **Extended Front Panels**  
Front panels are extended to 3/4 overlap coverage where LE protection is doubled
- Shoulder Padding**  
Premium foam padding in the shoulders provides extra comfort.
- Apron Pocket**
- Velcro Closure**
- Hanging Tabs on Skirt**

### ADDITIONAL CUSTOMIZATION OPTIONS

- Thyroid Collar
- Wings
- Extra Pocket
- Suspenders (skirt)
- Embroidery
- Logo
- Sleeve

| Core Material | Weight (0.25mm) | Thickness Options  |
|---------------|-----------------|--------------------|
| Lead          | 3400 gsm        | 0.25, 0.50mm       |
| Prolite       | 2850 gsm        | 0.25, 0.35, 0.50mm |



FOR SIZING AND COLOR OPTIONS SEE PAGES 55-60.



\*CE applies to Lead and Prolite in 0.25, 0.35, and 0.50mm. Prolite is limited to 110kVp diagnostic range.



FULL PROTECTION



# BALANCE™ FLEXBACK VEST AND SKIRT APRON

## APRON FEATURES

- +** **Adjustable Buckle Set**  
Adjustable buckle closures ensure a secure fit.
- +** **Flexback Belt**  
Wraps internally around your back and securely fastens in the front.
- +** **Non-slip Backing (skirt interior)**  
Reduces slippage and improves fit.
- +** **Seam-to-seam Protection**  
Seam-to-seam, standardized protection in Front (no change in protection at overlap)
- +** **Shoulder Padding**  
Premium foam padding in the shoulders provides extra comfort.
- +** **Apron Pocket**
- +** **Velcro Closure**
- +** **Hanging Tabs on Skirt**

### ADDITIONAL CUSTOMIZATION OPTIONS

- Thyroid Collar
- Wings
- Extra Pocket
- Suspenders (skirt)
- Embroidery
- Logo
- Sleeve

Core Material      Weight (0.25mm)      Thickness Options

|           |         |          |                    |
|-----------|---------|----------|--------------------|
| <b>Pb</b> | Lead    | 3400 gsm | 0.25, 0.50mm       |
| <b>P</b>  | Prolite | 2850 gsm | 0.25, 0.35, 0.50mm |



FOR SIZING AND COLOR OPTIONS SEE PAGES 55-60.




\*CE applies to Lead and Prolite in 0.25, 0.35, and 0.50mm. Prolite is limited to 110kVp diagnostic range.

MODEL: [AP-VFB]<sub>(vest only)</sub> [AP-SKIRT]<sub>(skirt only)</sub> [AP-VSFB]<sub>(vest and skirt combo)</sub>

The Vest Skirt Flexback shifts 70% of the weight directly to the hips with added support of an adjustable, elastic back relief belt. The Balance configuration features **seam-to-seam, single LE protection with standard overlapping panels** (see below for lead configuration). The skirt comes standard with an anti-slip fabric in the back (interior) to keep your skirt in place.

### Example of a 0.50mm LE Balance Vest Skirt Apron Lead Configuration



## DENTAL APRON



|                       |  |
|-----------------------|--|
| <b>MODEL</b>          | AP-DA  |
| <b>CORE MATERIALS</b> | Standard Lead or Prolite (see page 8)            |
| <b>THICKNESS</b>      | 0.50, 0.35, 0.25mm (standard) LE                 |
| <b>CLOSURE</b>        | None   |
| <b>SIZING</b>         | Adult: 28" W x 31" L<br>Pediatric: 24" W x 23" L |
| <b>CUSTOMIZATIONS</b> | Embroidery and Logo                              |
| <b>ADD-ONS</b>        | Thyroid Collar (included)                        |
| <b>COLORS</b>         | Refer to fabric colors on pages 55-60            |
| <b>CE</b>             | Yes  |

Features a non-slip microfiber material that ensures the apron lays in place over the patient. There are no closures. The thyroid collar is included.



## DENTAL PANORAMIC PONCHO



|                       |  |
|-----------------------|--|
| <b>MODEL</b>          | AP-DP  |
| <b>CORE MATERIALS</b> | Standard Lead or Prolite (see page 8)  |
| <b>THICKNESS</b>      | 0.50, 0.35, 0.25mm (standard) LE   |
| <b>CLOSURE</b>        | Velcro   |
| <b>SIZING</b>         | Adult: 24" W x 27" L x 16" L (front)<br>Pediatric: 18"W x 22"L (front) 12"L (back panel) |
| <b>CUSTOMIZATIONS</b> | Embroidery and Logo  |
| <b>ADD-ONS</b>        | Thyroid Collar   |
| <b>COLORS</b>         | Refer to fabric colors on pages 55-60  |
| <b>CE</b>             | Yes  |

Drapes over the body to provide front and back protection. Similar to the dental cape but the velcro closure is on the shoulder area. The Dental Panoramic Poncho comes in two sizes.



## ORIGINAL THYROID COLLAR



|                       |  |
|-----------------------|--|
| <b>MODEL</b>          | AP-TCO   |
| <b>CORE MATERIALS</b> | Standard Lead or Prolite (see page 8)              |
| <b>THICKNESS</b>      | 0.50 (standard), 0.35, 0.25mm LE                   |
| <b>CLOSURE</b>        | Velcro (shown), Magnets or Buckle                  |
| <b>SIZING</b>         | SM: 13-14"<br>MD: 15-16"<br>LG: 17"+<br>XL: Custom |
| <b>CUSTOMIZATIONS</b> | Embroidery and Logo                                |
| <b>ATTACHMENTS</b>    | Loose, Clip, Keyring, Tethered                     |
| <b>COLORS</b>         | Refer to fabric colors on pages 55-60              |
| <b>CE</b>             | Yes  |

The original style thyroid collar fits lower on the neck, offering coverage of the thyroid against scatter radiation. Velcro strips on the end allow for neck size adjustability. This streamlined style is our most popular.



## STANDARD THYROID COLLAR



|                       |  |
|-----------------------|--|
| <b>MODEL</b>          | AP-TCS   |
| <b>CORE MATERIALS</b> | Standard Lead or Prolite (see page 8)              |
| <b>THICKNESS</b>      | 0.50 (standard), 0.35, 0.25mm LE                   |
| <b>CLOSURE</b>        | Velcro, Magnets or Buckle (shown)                  |
| <b>SIZING</b>         | SM: 13-14"<br>MD: 15-16"<br>LG: 17"+<br>XL: Custom |
| <b>CUSTOMIZATIONS</b> | Embroidery and Logo                                |
| <b>ATTACHMENTS</b>    | Loose, Clip, Keyring, Tethered                     |
| <b>COLORS</b>         | Refer to fabric colors on pages 55-60              |
| <b>CE</b>             | Yes  |

The standard radiation thyroid collar is designed to fit higher along the neck while providing complete coverage of the thyroid against x-ray radiation.





## PROGUARD THYROID COLLAR



|                       |  |
|-----------------------|--|
| <b>MODEL</b>          | AP-TCG                                   |
| <b>CORE MATERIALS</b> | Standard Lead or Prolite (see page 8)    |
| <b>THICKNESS</b>      | 0.50 (standard), 0.35, 0.25mm LE         |
| <b>CLOSURE</b>        | Velcro                                   |
| <b>SIZING</b>         | SM: 21"<br>MD: 22"<br>LG: 23"<br>XL: 24" |
| <b>CUSTOMIZATIONS</b> | Embroidery and Logo                      |
| <b>ATTACHMENTS</b>    | Loose, Clip, Keyring, Tethered           |
| <b>COLORS</b>         | Refer to fabric colors on pages 55-60    |
| <b>CE</b>             | Yes                                      |

Featuring a form-fitting design and velcro closure, the Proguard thyroid collar provides front and back protection.



## MAMMOGRAPHY THYROID COLLAR



|                       |                                       |
|-----------------------|---------------------------------------|
| <b>MODEL</b>          | AP-TCM                                |
| <b>CORE MATERIALS</b> | Standard Lead or Prolite (see page 8) |
| <b>THICKNESS</b>      | 0.50 (standard), 0.35, 0.25mm LE      |
| <b>CLOSURE</b>        | Velcro (shown) or Buckle              |
| <b>SIZING</b>         | One Size: 29" L, 4" Bib               |
| <b>CUSTOMIZATIONS</b> | Embroidery and Logo                   |
| <b>ATTACHMENTS</b>    | Loose, Clip, Keyring, Tethered        |
| <b>COLORS</b>         | Refer to fabric colors on pages 55-60 |
| <b>CE</b>             | Yes                                   |

Don't let the standard or original thyroid collar get all the attention. The mammography's sleek profile provides a detachable alternative to "sewn-in" thyroid collars while allowing for adjustability. Additionally, the mammography thyroid collar serves as a valuable option for patient protection.



## THYROID COLLAR COVERS

|                       |                   |
|-----------------------|-------------------|
| <b>MODEL</b>          | SG-TCCD           |
| <b>CORE MATERIALS</b> | None              |
| <b>MATERIAL</b>       | SMS               |
| <b>SIZING</b>         | One size fits all |
| <b>FEATURES</b>       | Disposable        |
| <b>UOM QTY</b>        | 50 per box        |
|                       |                   |
|                       |                   |
|                       |                   |
|                       |                   |

Our disposable thyroid collar covers do not provide protection against x-ray, but are a useful way to encourage the use of thyroid collars by keeping them clean. The thyroid collar covers fit a wide variety of thyroid collar styles, will prevent damage from make-up and prevent the build-up of filth. Extend the life of your thyroid collars by covering them with ScatterGuard's hygienic, disposable thyroid collar cover.



## PROGUARD SMART CAP



|                       |                                       |
|-----------------------|---------------------------------------|
| <b>MODEL</b>          | AP-PGCAP                              |
| <b>CORE MATERIALS</b> | Standard Lead or Prolite (see page 8) |
| <b>THICKNESS</b>      | 0.50 (standard), 0.35, 0.25mm LE      |
| <b>CLOSURE</b>        | Beaded Adjustable Strap               |
| <b>SIZING</b>         | SM-M: 21-22"<br>L-XL: 23-25"          |
| <b>CUSTOMIZATIONS</b> | Embroidery and Logo                   |
| <b>COLORS</b>         | Refer to fabric colors on pages 55-60 |
| <b>CE</b>             | Yes                                   |
|                       |                                       |
|                       |                                       |

Need head protection that is breathable, light-weight, and comfortable? The Proguard Smart Cap features breathable CoolGuard spandex fabric on the top while providing protection around the sides of your head.



## PAGE CAP



|                       |                                       |
|-----------------------|---------------------------------------|
| <b>MODEL</b>          | AP-PAGE                               |
| <b>CORE MATERIALS</b> | Standard Lead or Prolite (see page 8) |
| <b>THICKNESS</b>      | 0.50 (standard), 0.35, 0.25mm LE      |
| <b>CLOSURE</b>        | Velcro                                |
| <b>SIZING</b>         | S: 21"<br>M: 22"<br>L: 23"<br>XL: 24" |
| <b>CUSTOMIZATIONS</b> | Embroidery and Logo                   |
| <b>PROTECTION</b>     | Side Protection                       |
| <b>COLORS</b>         | Refer to fabric colors on pages 55-60 |
| <b>CE</b>             | Yes                                   |

The Page Cap is designed to provide protection around the sides of your head, while the top contains no lead. Alternatively, similar to a scrub cap, the Head Cap provides x-ray protection around the entire head. Velcro tab closures allow for an adjustable fit on both styles. \*Sizes based on head circumference.



## HEAD CAP



|                       |                                       |
|-----------------------|---------------------------------------|
| <b>MODEL</b>          | AP-CAP                                |
| <b>CORE MATERIALS</b> | Standard Lead or Prolite (see page 8) |
| <b>THICKNESS</b>      | 0.50 (standard), 0.35, 0.25mm LE      |
| <b>CLOSURE</b>        | Velcro                                |
| <b>SIZING</b>         | S: 21"<br>M: 22"<br>L: 23"<br>XL: 24" |
| <b>CUSTOMIZATIONS</b> | None                                  |
| <b>PROTECTION</b>     | Top and Side Protection               |
| <b>COLORS</b>         | Refer to fabric colors on pages 55-60 |
| <b>CE</b>             | Yes                                   |

Similar to a scrub cap, the Head Cap provides x-ray protection around the entire head. Velcro tab closures allow for an adjustable fit on both styles. \*Sizes based on head circumference.





# SCOLIOSIS STOLE



|                       |  |
|-----------------------|--|
| <b>MODEL</b>          | AP-STOLE   |
| <b>CORE MATERIALS</b> | Standard Lead or Prolite (see page 8)                                |
| <b>THICKNESS</b>      | 0.50 (standard), 0.35, 0.25mm LE                                     |
| <b>CLOSURE</b>        | Buckle   |
| <b>SIZING</b>         | Pediatric: 15" x 4.75"<br>Medium: 18.5" x 6.75"<br>Large: 23.5" x 8" |
| <b>CUSTOMIZATIONS</b> | Embroidery and Logo  |
| <b>COLORS</b>         | Refer to fabric colors on pages 55-60                                |
| <b>CE</b>             | Yes  |

The scoliosis stole is primarily used for spinal imaging. The shield features velcro panels and an adjustable belt for proper positioning. The neck area does not contain lead, allowing complete visibility of the spine during x-ray, while offering protection to the breast area. This shield is great for scoliosis studies.



## GONAD/OVARIAN SHIELD



|                       |  |
|-----------------------|--|
| <b>MODEL</b>          | AP-GONAD/GONAD-SET   |
| <b>CORE MATERIALS</b> | Standard Lead or Prolite (see page 8)  |
| <b>THICKNESS</b>      | 0.50 (standard), 0.35, 0.25mm LE   |
| <b>CLOSURE</b>        | Velcro (standard) or Buckle  |
| <b>SIZING</b>         | Small - 5.75" x 4.5"<br>Medium - 7.6" x 6.4"<br>Large - 9.5" x 7.8"<br>Pediatric sizing also available (page 54) |
| <b>CUSTOMIZATIONS</b> | Embroidery and Logo  |
| <b>SPECIAL</b>        | Available as a set of 3 or individually  |
| <b>COLORS</b>         | Refer to fabric colors on pages 55-60  |
| <b>CE</b>             | Yes  |

This unique gonad/ovarian shield belt system comes complete with four interchangeable sizes which provide coverage from x-ray for a variety of needs. The adjustable Velcro or buckle also allows the user to fit almost any waist size. The gonad/ovarian belt system comes as a set of 3 sizes (adult or pediatric sets). They can also be purchased individually.



## APRON SLEEVE




|                       |                                       |
|-----------------------|---------------------------------------|
| <b>MODEL</b>          | AP-SLV/SLV2                           |
| <b>CORE MATERIALS</b> | Standard Lead or Prolite (see page 8) |
| <b>THICKNESS</b>      | 0.50, 0.35, 0.25mm (standard) LE      |
| <b>CLOSURE</b>        | Snap-on or Sewn-in                    |
| <b>SIZING</b>         | Small, Medium, Large, XL or Custom    |
| <b>CUSTOMIZATIONS</b> | Embroidery and Logo                   |
| <b>COLORS</b>         | Refer to fabric colors on pages 55-60 |
| <b>CE</b>             | Yes                                   |

This optional sleeve is available for many apron styles and protects the upper arm and shoulder from direct beam and scatter radiation. It can be attached as a snap-on or directly sewn-on. Please note that the apron sleeve is not fully attached at the armpit. An optional flap is available for additional protection/coverage (see customizations page for more info).

Sewn-in



Snap-on



## SHIN GUARDS



|                       |   |
|-----------------------|---|
| <b>MODEL</b>          | AP-SG   |
| <b>CORE MATERIALS</b> | Standard Lead or Prolite (see page 8)                         |
| <b>THICKNESS</b>      | 0.50 (standard), 0.35, 0.25mm LE                              |
| <b>CLOSURE</b>        | Velcro  |
| <b>SIZING</b>         | SM: 13" x 9"<br>M: 14" x 10"<br>L: 16" x 12"<br>XL: 18" x 24" |
| <b>CUSTOMIZATIONS</b> | Embroidery and Logo   |
| <b>COLORS</b>         | Refer to fabric colors on pages 55-60                         |
| <b>CE</b>             | Yes   |

When being exposed to ionizing radiation, all forms of radiation protection are necessary. Shin Guards protect against beam scatter when no table drapes are present. Optional foot covers may be added (AP-FOOT).



(AP-FOOT)

## HALF APRON



|                       |  |
|-----------------------|--|
| <b>MODEL</b>          | AP-HA  |
| <b>CORE MATERIALS</b> | Standard Lead or Prolite (see page 8)  |
| <b>THICKNESS</b>      | 0.50 (standard), 0.35, 0.25mm LE   |
| <b>CLOSURE</b>        | Velcro   |
| <b>SIZING</b>         | SM (pediatric): 10" x 10"<br>MD: 18" x 16"<br>LG: 24" x 20"<br>XL: 24" x 24" |
| <b>CUSTOMIZATIONS</b> | Embroidery and Logo  |
| <b>COLORS</b>         | Refer to fabric colors on pages 55-60  |

The half apron can provide extra radiation protection. Selling it separate also allows you to mix and match a smaller bottom half with a larger top piece or vice versa.





## CORSET BELT

|                 |   |
|-----------------|---|
| MODEL           | AP-CORSET                                 |
| PROTECTION      | None                                      |
| SIZING (FEMALE) | Small - 4XL (Use apron waist measurement) |
| SIZING (MALE)   | Small - 4XL (Use apron waist measurement) |
| COLORS          | Black                                     |
|                 |   |
|                 |   |
|                 |   |
|                 |   |
|                 |   |
|                 |   |

Our corset belt is a great apron accessory for added back support when wearing a radiation protective apron. The finger pull tab and pulley system along the spine provides for simple adjustment and even tension.

**\*The corset belt is worn under the apron.**  
**\*\* Use apron waist measurements on page 54 to determine your size.**



## BACK RELIEF BELT

|                 |                                       |
|-----------------|---------------------------------------|
| MODEL           | AP-BLT                                |
| PROTECTION      | None                                  |
| CLOSURE         | Buckle                                |
| SIZING (FEMALE) | XS: 24-26"   S: 26-28"                |
|                 | M: 28-31"   L: 31-33"                 |
|                 | XL: 34-38"                            |
| SIZING (FEMALE) | XS: 30-33"   S: 33-36"                |
|                 | M: 37-40"   L: 41-44"                 |
|                 | XL: 44-48"                            |
| CUSTOMIZATIONS  | None                                  |
| COLORS          | Refer to fabric colors on pages 55-60 |

The back relief belt is constructed of neoprene bonded elastic. It features high quality components that do not lose elasticity. The belt also features a non-slip adjustable buckle to ensure a secure fit. The back relief belt attaches separately from the apron.



# PEDIATRIC PROTECTION



## PEDIATRIC TIE APRON



|                       |  |
|-----------------------|--|
| <b>MODEL</b>          | AP-TA-CH   |
| <b>CORE MATERIALS</b> | Lead or Prolite (see page 8)   |
| <b>THICKNESS</b>      | 0.50 (standard), 0.35, 0.25mm LE   |
| <b>CLOSURE</b>        | Tie  |
| <b>SIZING</b>         | Small: 12 x 24" (Ages 1-3)<br>Medium: 14 x 26" (Ages 4-6)<br>Large: 16 x 28" (Ages 7-10) |
| <b>CUSTOMIZATIONS</b> | Embroidery and Logo  |
| <b>COLORS</b>         | Refer to fabric colors on pages 55-60  |
| <b>CE</b>             | Yes  |

The Tie Apron is an easy-to-use, affordable option for x-ray protection. Utilizing a tie closure, the Tie Apron's criss crossing straps ensure a secure and comfortable fit. Premium foam padding in the shoulders provides extra comfort.



## PEDIATRIC VELCRO ADJUSTABLE APRON



|                       |  |
|-----------------------|--|
| <b>MODEL</b>          | AP-VA-CH   |
| <b>CORE MATERIALS</b> | Standard Lead or Prolite (see page 8)  |
| <b>THICKNESS</b>      | 0.50 (standard), 0.35, 0.25mm LE   |
| <b>CLOSURE</b>        | Velcro   |
| <b>SIZING</b>         | Small: 12 x 24" (Ages 1-3)<br>Medium: 14 x 26" (Ages 4-6)<br>Large: 16 x 28" (Ages 7-10) |
| <b>CUSTOMIZATIONS</b> | Embroidery and Logo  |
| <b>COLORS</b>         | Refer to fabric colors on pages 55-60  |
| <b>CE</b>             | Yes  |

The Velcro Adjustable apron has crisscrossing back straps that provide back support while securing the apron around your waist. Wide back straps prevent the apron from slipping around your shoulder and back area.





## PEDIATRIC HALF APRON



|                       |                                       |
|-----------------------|---------------------------------------|
| <b>MODEL</b>          | AP-HA-CH                              |
| <b>CORE MATERIALS</b> | Standard Lead or Prolite (see page 8) |
| <b>THICKNESS</b>      | 0.50 (standard), 0.35, 0.25mm LE      |
| <b>CLOSURE</b>        | Velcro                                |
| <b>SIZING</b>         | 10" x 10"                             |
| <b>CUSTOMIZATIONS</b> | Embroidery and Logo                   |
| <b>COLORS</b>         | Refer to fabric colors on pages 55-60 |
| <b>CE</b>             | Yes                                   |

The half apron can provide extra radiation protection. Selling it separate also allows you to mix and match a smaller bottom half with a larger top piece or vice versa.



## PEDIATRIC DENTAL APRON



|                       |                                       |
|-----------------------|---------------------------------------|
| <b>MODEL</b>          | AP-DA-CH                              |
| <b>CORE MATERIALS</b> | Standard Lead or Prolite (see page 8) |
| <b>THICKNESS</b>      | 0.50, 0.35, 0.25mm (standard) LE      |
| <b>CLOSURE</b>        | None                                  |
| <b>SIZING</b>         | Pediatric: 24" W x 23" L              |
| <b>CUSTOMIZATIONS</b> | Embroidery and Logo                   |
| <b>SPECIAL</b>        | Thyroid Collar (included)             |
| <b>COLORS</b>         | Refer to fabric colors on pages 55-60 |
| <b>CE</b>             | Yes                                   |

Features a non-slip microfiber material that ensures the apron lays in place over the patient. There are no closures. The thyroid collar is included.



## PEDIATRIC THYROID COLLAR



|                       |  |
|-----------------------|--|
| <b>MODEL</b>          | AP-TCP   |
| <b>CORE MATERIALS</b> | Standard Lead or Prolite (see page 8)                |
| <b>THICKNESS</b>      | 0.50 (standard), 0.35, 0.25mm LE                     |
| <b>CLOSURE</b>        | Velcro   |
| <b>SIZING</b>         | Infant: 10" (Ages 2-4)<br>Pediatric: 12" (Ages 5-10) |
| <b>CUSTOMIZATIONS</b> | Embroidery and Logo                                  |
| <b>ATTACHMENTS</b>    | Loose, Clip, Keyring, Tethered                       |
| <b>COLORS</b>         | Refer to fabric colors on pages 55-60                |
| <b>CE</b>             | Yes  |

Available in either the standard or original style thyroid collars, our Pediatric Thyroid collar is designed to protect children ages 2 to 10 from harmful x-ray.



## PEDIATRIC GONAD SHIELD SET



|                       |   |
|-----------------------|---|
| <b>MODEL</b>          | AP-GONADSET2  |
| <b>CORE MATERIALS</b> | Standard Lead or Prolite (see page 8)                                   |
| <b>THICKNESS</b>      | 0.50 (standard), 0.35, 0.25mm LE  |
| <b>CLOSURE</b>        | Velcro (standard) or Buckle   |
| <b>SIZING</b>         | 1 Pediatric: 3.75" x 3.75"<br>1 Child: 5" x 4"<br>1 Small: 5.75" x 4.5" |
| <b>CUSTOMIZATIONS</b> | Embroidery and Logo   |
| <b>COLORS</b>         | Refer to fabric colors on pages 55-60                                   |
| <b>CE</b>             | Yes   |

This unique pediatric gonad/ovarian shield belt system comes complete with three interchangeable sizes which provide coverage from x-ray for a variety of needs. The adjustable Velcro or optional buckle (shown) also allows the user to fit almost any waist size. The gonad/ovarian belt system comes as a set of 3 sizes, or can be purchased individually.



## PEDIATRIC SCOLIOSIS STOLE



|                       |                                       |
|-----------------------|---------------------------------------|
| <b>MODEL</b>          | AP-STOLE-CH                           |
| <b>CORE MATERIALS</b> | Standard Lead or Prolite (see page 8) |
| <b>THICKNESS</b>      | 0.50 (standard), 0.35, 0.25mm LE      |
| <b>CLOSURE</b>        | Buckle                                |
| <b>SIZING</b>         | Pediatric: 15" x 4.75"                |
| <b>CUSTOMIZATIONS</b> | Embroidery and Logo                   |
| <b>COLORS</b>         | Refer to fabric colors on pages 55-60 |
| <b>CE</b>             | Yes                                   |

The scoliosis stole is primarily used for spinal imaging. The shield features velcro panels and an adjustable belt for proper positioning. The neck area does not contain lead, allowing complete visibility of the spine during x-ray, while offering protection to the breast area. This shield is great for scoliosis studies.



## TUGGA BABY



|                     |                                     |
|---------------------|-------------------------------------|
| <b>MODEL</b>        | RE-BABY                             |
| <b>FRAME TYPE</b>   | Goggles                             |
| <b>FRONT LENS</b>   | 0.75mm Pb Glass                     |
| <b>WEIGHT</b>       | 42g                                 |
| <b>FRAME SIZE</b>   | 46-18                               |
| <b>MATERIAL</b>     | N/A                                 |
| <b>SIDE SHIELDS</b> | N/A                                 |
| <b>RX OPTIONS</b>   | N/A                                 |
| <b>MIN PRG PD</b>   | N/A                                 |
| <b>CE</b>           | Yes                                 |
| <b>COLORS</b>       | Shiny Pink, Shiny Blue, Shiny Green |

Tugga Baby pediatric glasses provide essential protection during medical fluoroscopy and feature a universal fit with its adjustable headband.





# APRON STORAGE

## EXTEND THE LIFE OF YOUR APRON


When not in use, garments should be placed on a hanger or stored using Proguard's apron storage equipment. Avoid contact with hot surfaces or exposure to direct sunlight. Also avoid contact with sharp objects, as these may damage the garment's protective core material. For transportation, garments should be carefully rolled up, not folded, and placed into a box or Proguard duffle bag.

### WALL RACK

MODEL: [ AS-WALL-L/R ]

CAPACITY: 5 APRONS PER CONSOLE

Prolong the life of your radiation protective garment with our apron wall rack console. Each console holds up to 5 protective aprons with the combined left and right consoles storing up to 10 aprons. The wall console features hooks on the top of the hangers for lead skirts. \*Does not include mounting hardware.




### PEG RACK

MODEL: [ AS-PEG ]

CAPACITY: 1-10 APRONS

The apron Peg rack is made of steel construction and easily screws into wall studs. It can store up to 10 aprons but comes in a variety of sizes that accommodate a range of lead aprons. This affordable wall storage for aprons is essential to prolong the life of your protective aprons. \*Mounting hardware is not included.



| Model  | # of Pegs | # of Aprons | Dimensions |
|--------|-----------|-------------|------------|
| PEG-2  | 2         | 1           | 5 x 3"     |
| PEG-4  | 4         | 2           | 10 x 3"    |
| PEG-6  | 6         | 3           | 15 x 3"    |
| PEG-8  | 8         | 4           | 20 x 3"    |
| PEG-10 | 10        | 5           | 25 x 3"    |
| PEG-14 | 14        | 7           | 35 x 3"    |
| PEG-18 | 18        | 9           | 45 x 3"    |
| PEG-20 | 20        | 10          | 50 x 3"    |

### PORTA RACK

MODEL: [ AS-PORTRK ]

CAPACITY: 8 APRONS

Ideal for use in the O.R., Specials, and Cath Labs, this portable apron storage unit is great quality and affordable. The 24" wheel base with double casters allows this completely loaded rack to glide smoothly while being moved. The porta rack features moving hanger arms that swing forward and backwards for easy apron selection.



## MIGHTY MAX RACK

MODEL: [ AS-MAX-MOB ]

CAPACITY: 10 APRONS

The Mighty-Max mobile apron rack's solid steel design holds up to 10 aprons. The pedestal base has rolling wheels for easy portability. Convenience hooks on each hanger can hold lead skirts.



## APRON HANGER

MODEL: [ AS-HANGER ]

CAPACITY: 40LBS

Constructed out of aluminum for durability and longevity this apron hanger is used for the storage of radiation protection aprons. It features hook style ends that are designed to hold lead skirts.



## VALET RACK

MODEL: [ AS-VALET/VALET-VS ]

CAPACITY: UP TO 15 APRONS/10 APRONS(VS)

Industrial quality steel construction and easy roll casters make these mobile apron racks an excellent portable storage option. The rack provides storage capacity for up to 15 aprons! Casters are heavy-duty, medical-grade (AS-Valet-VS)



## APRON DUFFLE BAG

MODEL: [ AS-DUFFLE ]

CAPACITY: UP TO 2 APRONS


Our duffle bag protects your aprons while traveling between facilities or out of the country. Protect and store your aprons for transport now with our affordable, high quality duffle bag. Features: 10" Velcro Opening, Carrying Handle & Shoulder Strap.

For all color options, refer to pages 55-60.



# APRON MEASUREMENTS & SIZING

## CHEST MEASUREMENT




Measure around the largest circumference of the bust area, keeping the tape parallel to the floor with the client's arms at their sides.

## WAIST MEASUREMENT



Measure around the smallest circumference of the torso. If you measure at waist level it may be inaccurate.

## HIP MEASUREMENT




Measure around the largest circumference of the hips and buttocks with the tape parallel to the floor: generally 8-9 inches down from the waist.

## FRONT APRON LENGTH MEASUREMENT




Measure from the top of the shoulder over the bust to the desired length. Do not contour the tape measure in under the bust.

## VEST LENGTH MEASUREMENT



Measure from the top of the shoulder over the bust to the desired length. Do not contour the tape measure under the bust.

## SKIRT LENGTH MEASUREMENT




Measure from the top of the pant line (waist) and over the outside of the hip down to the desired length.



**ALL FRONTAL & WRAP APRONS**

**TWO-PIECE VEST & SKIRT APRONS**



**MEN**  
  
**WOMEN**  
  
**KIDS\***

|              | XS   | Small | Med  | Large | XL    | 2X   | 3X   |
|--------------|------|-------|------|-------|-------|------|------|
| Height       | 5'2" | 5'4"  | 5'8" | 5'11" | 6'2"+ |      |      |
| Chest        | 39"  | 42"   | 45"  | 48"   | 53"   | 57"  | 60"  |
| Waist        | 32   | 34    | 38   | 42    | 46    | 50   | 54   |
| Hips         | 34   | 36    | 40   | 44    | 48    | 52   | 56   |
| Front Length | 34   | 34    | 37   | 40.5  | 42.5  | 42.5 | 45.5 |
| Front Width  | 19   | 21    | 22   | 24    | 26    | 28   | 30   |
| Vest Length  | 22   | 22    | 24   | 25    | 25    | 25   | 28   |
| Skirt Length | 20   | 20    | 20   | 21    | 21    | 22   | 22   |
| WABR Length  | 36   | 36    | 38   | 40.5  | 42.5  | 42.5 | 45.5 |

|              | XS   | Small | Med  | Large | XL    | 2X    | 3X   |
|--------------|------|-------|------|-------|-------|-------|------|
| Height       | 5'1" | 5'2"  | 5'5" | 5'8"  | 5'10" | 5'11" |      |
| Chest        | 36"  | 38"   | 40"  | 43"   | 46"   | 49"   | 52"  |
| Waist        | 25   | 27    | 29   | 32    | 35    | 40    | 44   |
| Hips         | 35   | 37    | 39   | 42    | 46    | 48    | 52   |
| Front Length | 33   | 34    | 36   | 38.5  | 40.5  | 42.5  | 42.5 |
| Front Width  | 18   | 19    | 20   | 21    | 23    | 24    | 26   |
| Vest Length  | 19   | 22    | 22   | 22    | 22    | 23    | 23   |
| Skirt Length | 20   | 20    | 22   | 22.5  | 22.5  | 22.5  | 24   |
| WABR Length  | 33   | 34    | 37   | 40.5  | 42.5  | 42.5  | 44.5 |

|            | Small    | Medium   | Large     |
|------------|----------|----------|-----------|
| Dimensions | 12 x 24" | 14 x 26" | 16 x 28"  |
| Ages       | Ages 1-3 | Ages 4-6 | Ages 7-10 |


\*Velcro Adjustable & Tie Aprons

# FABRIC COLORS

## LIMITED EDITION PRINTS\*



Ankara Cloth  
#2000



Blue Marble  
#2010




Body Positivity  
#2015



Exhaust-ed  
#2020




Galactic Kitties  
#2025



Geo Tribe  
#2030




Grind my Gears  
#2035



Hex Yeah!  
#2040




Medically Necessary  
#2045



Midnight Jungle  
#2050




Party 'Till You're Dead  
#2070



Reel Time  
#2075



Stable Genius  
#2080



Stardate: 2022  
#2022



Starry Night  
#2090



Supernova  
#2095



This Lead is Bananas  
#2005




Tropical Smoothie  
#2100






## COOLGUARD PERFORMANCE FABRICS\*




**Coolguard Air**  
#900

A cushiony mesh material that enables natural air ventilation. Interior only.



**Coolguard Spandex**  
#915

Sweat wicking, breathable microfiber with quick dry technology.




**Coolguard Mesh**  
#910

Increases air flow on back closure of apron. For frontal apron back closures only.



**Indigo**  
#925

A soft, PVC nylon that is fluid-proof, stain resistant, bacteria & fungus resistant. Our lightest-weight fabric!



**"Black Silk"**  
#920

Ultra-soft & breathable polyester that is fluid proof & stain resistant.




**Anti-slip fabric**  
#950

A high-quality, textured fabric placed on the inside of your apron to prevent slippage.

## SATINGUARD ANTI-MICROBIAL FABRICS\*




Featuring stain resistant and anti-microbial technology only available from Proguard. Our soft poly material will resist a wide variety of bacteria for several years.




Matte Black  
#800



Matte Bordeaux Red  
#824



Matte Cyan  
#805



Matte Light Gray  
#815



Matte Dark Gray  
#813



Matte Lilac  
#817



Matte Navy  
#820



Matte Fuschia  
#810



Matte Red  
#825



Matte Royal  
#802



Matte Dark Teal  
#828



Matte Light Teal  
#830

\*Upcharge applies



# NYLONS

Unless specified otherwise, nylons come standard 210 denier.



Black  
#100



Burgundy  
#105



Forest Green  
#110




Gray  
#115




Green  
#120




Hot Pink  
#125




Marine Blue  
#135



Navy  
#140




Neon Green  
#145



Orange  
#150




Purple  
#155



Raspberry  
#160



Red  
#165



Rose  
#170




Royal  
#175



Tan  
#180



Teal  
#185



Yellow  
#190

# BINDING COLORS



Purple



Burgundy



Red



Orange



Yellow



Neon Green



Green



Forest Green



Teal



Marine Blue



Royal



Navy Blue



Black



Gray



Tan



Rose




Hot Pink




Raspberry


## RIPSTOP




Black  
#200




Brown  
#205




Burgundy  
#210




Chocolate  
#215




Forest Green  
#220




Green  
#225




Navy  
#230



Purple  
#235



Red  
#240




Royal Blue  
#245




Silver  
#250

## WIPE-AWAY VINYL\*



Black  
#300



Silver  
#320



Burgundy  
Starlight  
#305

\*Upcharge applies

# CLASSIC PRINTS



70's Flowers  
#400



Animal Farm  
#405



Argyle Green  
#410




Asian Garden  
#415




Black Lace  
#417




Blue Animal River  
#420



Blue Flame  
#425



Blue Tranquility  
#430




Brown Animal Spots  
#435




Butterfly  
#440



Camo - Backwoods  
#443



Camo - Cosmo  
#450



Camo - Digital Grey  
#455



Camo - Pink  
#465




Camo - Small Blues  
#470



Camo - Small Pinks  
#475



Camo - Yellow  
#480



Chicklets Blue  
#490



Circles  
#495



Crayon Hearts  
#500



Divebar  
#502



Dog Paw  
#505



Fire Tongues  
#510



Fleur De Lis  
#512



Fruity Pebbles  
#515




Galaxy - Multi  
#520



Galaxy - Royal  
#525



Giraffe - Brown  
#530



Green Neon Tie-Dye  
#540



Hawaiian  
#545



Hawaii Postcard  
#550



Leopard Print  
#560




Lu Lu Grape  
#565



Mamba - Blue  
#575



Mamba - Pink  
#580



Marble Splash  
#585



# CLASSIC PRINTS




Mardi Gras  
#590



Mystery - Multi  
#595



Painted Dragon  
#600



Paint Splash  
#605



Paw Prints  
#610



Plaid - Blue  
#615



Plaid - Green  
#620



Plaid - Light Pink  
#625



Plaid - Purple  
#630



Plaid - Tan  
#635



Plaid - Teal  
#640




Purple Neon  
Tie-Dye #645



Rainbow  
#650



Red Dalmatian  
#655



Sheep  
#665




Skulls - Gray  
#670



Skulls - Red  
#675



Snake Skin  
#680



Star Gazer  
#685




Stars  
#690




Steel  
#695




Surf Wax  
#705



Teal Splash  
#710



Tiger  
#715



Zebra  
#720

\*Upcharge applies









## CONTACT US

---

**Telephone: +44 (0)1606 331119**

**Email: [info@apccv.com](mailto:info@apccv.com)**

**Website: [www.xrayaprons.co.uk](http://www.xrayaprons.co.uk)**

APC Cardiovascular Ltd.  
Whitehall Business Centre  
75 School Lane, Hartford CW8 1PF

Registered in England No. 2423606.