APRON MEASUREMENTS

CHEST MEASUREMENT



Measure around the largest circumference of the bust area, keeping the tape parallel to the floor with the client's arms at their sides.

WAIST MEASUREMENT



Measure around the smallest circumference of the torso. If you measure at waist level it may be inaccurate.

HIP MEASUREMENT



Measure around the largest circumference of the hips and buttocks with the tape parallel to the floor: generally 8-9 inches down from the waist.

FRONT LENGTH MEASUREMENT



Measure from the top of the shoulder over the bust to the desired length. Do not contour the tape measure in under the bust.

VEST LENGTH MEASUREMENT




Measure from the top of the shoulder over the bust to the desired length. Do not contour the tape measure under the bust.

SKIRT LENGTH MEASUREMENT



Measure from the top of the pant line (waist) and over the outside of the hip down to the desired length.

| | | | | A | P | | 3 | | |
|---|------------------|----------------|------------|------------|------------|-----------------|----------------|----------------|-----------------|
| | | Height | Chest | Waist | Hips | Front Length | Front Width | Vest Length | Skirt Length |
| | Small Medium | 5'4" 5'8" | 42" 45" | 34" 38" | 36" 40" | 34" 37" | 21" 22" | 22" 24" | 20" 20" |
| 5 | Large X-Large | 5'11" 6'2"+ | 48" 53" | 42" 46" | 44" 48" | 40.5" 42.5" | 24" 26" | 25" 25" | 21" 21" |
| | | | | | | P | | | |



PRETECH

APRON SIZE CHART (UNISEX)

Use the measurements from our apron size chart to determine your closest apron size. Need help finding your size? Contact us at **561.627.9769** or email **sales@protechmed.com** for additional assistance.

Apollo Quickship aprons are only available in standard sizing (no custom, made-to-measure).



| | XS | Small | Med | Large | XL | 2X | 3X |
|--------------|------|-------|------|-------|-------|------|------|
| Height | 5'2" | 5'4" | 5'8″ | 5'11" | 6'2"+ | | |
| Chest | 39" | 42" | 45" | 48" | 53″ | 57" | 60" |
| Waist | 32 | 34 | 38 | 42 | 46 | 50 | 54 |
| Hips | 34 | 36 | 40 | 44 | 48 | 52 | 56 |
| Front Length | 34 | 34 | 37 | 40.5 | 42.5 | 42.5 | 45.5 |
| Front Width | 19 | 21 | 22 | 24 | 26 | 28 | 30 |
| Vest Length | 22 | 22 | 24 | 25 | 25 | 25 | 28 |
| Skirt Length | 20 | 20 | 20 | 21 | 21 | 22 | 22 |
| WABR Length | 36 | 36 | 38 | 40.5 | 42.5 | 42.5 | 45.5 |

ALL FRONTAL & TWO-PIECE VEST **& SKIRT APRONS** WRAP APRONS Front Width Vest Length Chest Highest point of shoulder by neck, **Front Length** to waist. Add 2" for Measure from overlap highest point of shoulder to the knee or desired length Waist 3" Overlap Hips **Skirt Length** Measure from waistline to the knee or desired length

APRON SIZE CHART (MALE)

Use the measurements from our apron size chart to determine your closest apron size. Need help finding your size? Contact us at **561.627.9769** or email **sales@protechmed.com** for additional assistance.

Looking for a made-to-measure apron? Contact us and provide your measurements for a fully custom apron!



| | | XS | Small | Med | Large | XL | 2X | 3X | |
|-----|--------------|------|-------|------|-------|-------|-------|------|--|
| | Height | 5'1" | 5'2" | 5'5″ | 5'8″ | 5'10" | 5′11″ | | |
| M | Chest | 36" | 38" | 40" | 43" | 46″ | 49" | 52" | |
| 6 | Waist | 25 | 27 | 29 | 32 | 35 | 40 | 44 | |
| | Hips | 35 | 37 | 39 | 42 | 46 | 48 | 52 | |
| | Front Length | 33 | 34 | 36 | 38.5 | 40.5 | 42.5 | 42.5 | |
| - 1 | Front Width | 18 | 19 | 20 | 21 | 23 | 24 | 26 | |
| | Vest Length | 19 | 22 | 22 | 22 | 22 | 23 | 23 | |
| 1 | Skirt Length | 20 | 20 | 22 | 22.5 | 22.5 | 22.5 | 24 | |
| 11 | WABR Length | 33 | 34 | 37 | 40.5 | 42.5 | 42.5 | 44.5 | |

APRON SIZE CHART (FEMALE)

Use the measurements from our apron size chart to determine your closest apron size. Need help finding your size? Contact us at **561.627.9769** or email **sales@protechmed.com** for additional assistance.

Looking for a made-to-measure apron? Contact us and provide your measurements for a fully custom apron!



