



PROTECH LEADED EYEWEAR INC,
DBA PROTECH MEDICAL

1360 North Killian Drive, Unit 2
Lake Park, FL 33403

Email(USA): sales@protechmed.com
Email(Intl): global@protechmed.com
Phone: 561.627.9769
Web: www.protechmed.com

LEAD GARMENT WEIGHTS

FRONT PROTECTION APRONS

A. Frontal Apron Styles

Buckle Apron

Tie Apron

Back Relief Apron

Flexback Apron

Surgical Drop Apron

Velcro Adjustable Apron

Urology Apron



B. Average Frontal Apron Dimensions and Weights

Protech makes every effort to accurately record its apron weights. Because each category is averaged, the numbers we provide in the weight chart are subject to a +/- 15% weight variance. Minor length or width adjustments, including heavier fabrics, are not differentiated and instead, averaged together. Please consider this when referencing the weight charts.

AVERAGE FRONTAL 0.50MM LE APRON DIMENSIONS (IN) AND WEIGHTS (LBS)

SIZING	DIMENSIONS		WEIGHTS			
	LENGTH	WIDTH	LEAD	PROLITE	PREMIER	NANO
X-Small	34.75"	20.25"	7.17	5.95	5.59	4.87
Small	35.75"	21.00"	7.68	6.37	5.99	5.22
Medium	38.25"	22.50"	8.63	7.16	6.73	5.87
Large	41.50"	24.50"	10.35	8.59	8.07	7.04
X-Large	44.00"	27.00"	11.74	9.74	9.16	7.98
2X-Large	44.00"	28.25"	12.28	10.20	9.58	8.35
3X-Large	46.00"	31.25"	14.21	11.79	11.08	9.66



PROTECH LEADED EYEWEAR INC,
DBA PROTECH MEDICAL

1360 North Killian Drive, Unit 2
Lake Park, FL 33403

Email(USA): sales@protechmed.com
Email(Intl): global@protechmed.com
Phone: 561.627.9769
Web: www.protechmed.com

LEAD GARMENT WEIGHTS

FULL PROTECTION WRAPAROUND APRONS

A. Wraparound Apron Styles

Essential Wraparound
Back Relief Apron



Balance Wraparound
Back Relief Apron



Wraparound
Flexback Apron



Wraparound Flexback
Maternity Apron



B. Average Wraparound Apron Dimensions and Weights

Protech makes every effort to accurately record its apron weights. Because each category is averaged, the numbers we provide in the weight chart are subject to a +/- 15% weight variance. Minor length or width adjustments, including heavier fabrics, are not differentiated and instead, averaged together. Please consider this when referencing the weight charts.

AVERAGE WRAPAROUND 0.50MM LE APRON DIMENSIONS (IN) AND WEIGHTS (LBS)

SIZING	DIMENSIONS		WEIGHTS			
	LENGTH	WIDTH	LEAD	PROLITE	PREMIER	NANO
X-Small	-	-	12.44	10.33	9.70	8.46
Small	35.75"	19.25"	14.40	11.95	11.23	9.79
Medium	39.00"	19.50"	15.95	13.23	12.44	10.84
Large	41.50"	23.00"	17.52	14.54	13.67	11.91
X-Large	43.00"	24.00"	19.08	15.84	14.89	12.98
2X-Large	44.50"	27.00"	22.00	18.26	17.16	14.96
3X-Large	-	-	-	-	-	-



PROTECH LEADED EYEWEAR INC,
DBA PROTECH MEDICAL

1360 North Killian Drive, Unit 2
Lake Park, FL 33403

Email(USA): sales@protechmed.com
Email(Intl): global@protechmed.com
Phone: 561.627.9769
Web: www.protechmed.com

LEAD GARMENT WEIGHTS

FULL PROTECTION VEST SKIRT APRONS

A. Vest Skirt Apron Styles

Essential Vest Skirt Apron



Balance Vest Skirt Apron



Balance Skirt Flexback Apron



Essential Skirt Flexback Apron



Tri-tab Vest Skirt Apron



B. Average Vest Skirt Apron Dimensions and Weights

Protech makes every effort to accurately record its apron weights. Because each category is averaged, the numbers we provide in the weight chart are subject to a +/- 15% weight variance. Minor length or width adjustments, including heavier fabrics, are not differentiated and instead, averaged together. Please consider this when referencing the weight charts.

AVERAGE VEST SKIRT 0.50MM LE APRON DIMENSIONS (IN) AND WEIGHTS (LBS)

SIZING	DIMENSIONS				WEIGHTS			
	VEST LENGTH	VEST WIDTH	SKIRT LENGTH	SKIRT WIDTH	LEAD	PROLITE	PREMIER	NANO
X-Small	22"	19.5"	20"	17"	11.16	9.32	8.54	7.53
Small	22"	21.0"	20"	18"	12.74	10.64	9.75	8.60
Medium	24"	22.5"	20"	20"	14.75	12.32	11.29	9.96
Large	25"	24.0"	21"	22"	17.96	15.00	13.74	12.12
X-Large	25"	26.5"	21"	24"	19.56	16.33	14.96	13.20
2X-Large	25"	28.5"	21"	26"	21.46	17.92	16.42	14.48
3X-Large	28"	30.0"	21"	28"	23.84	19.91	18.24	16.09