


# PROGUARD **BALANCE**

BALANCE FLEXBACK VEST AND SKIRT APRON  
- 0.50MM LE LEAD CONFIGURATION -



## **BALANCE**


- **Standard Overlapping Panel**
- **Seam-to-seam, single LE protection**
- **Standard 0.25mm LE in back of skirt**
- **Custom and Standard Sizing available**

Dark gray areas indicate where panel(s) are placed to create desired protection. In this example, there is uniform 0.50mm LE protection from seam to seam (no protection changes at overlap). This is achieved through the use of additional, internal layers of lead.

**Example**  
**Order Configuration - 0.50/0.25mm LE**

# PROGUARD ESSENTIAL

ESSENTIAL FLEXBACK VEST AND SKIRT APRON  
- 0.50MM LE LEAD CONFIGURATION -



## ESSENTIAL

- 3/4 Overlapping Panel
- Desired protection at overlap where most sensitive organs are - helps with weight reduction
- Standard 0.25mm LE in back of skirt
- Standard Sizing Only

Dark gray areas indicate where vest panels overlap to create desired protection. In this example, two 0.25mm vest panels overlap to create 0.50mm LE protection over the body's most sensitive organs.

**Example**  
**Order Configuration - 0.50/0.25mm LE**