


PROGUARD **BALANCE**

BALANCE WRAPAROUND BACK RELIEF APRON
- 0.50MM LE LEAD CONFIGURATION -



BALANCE


- Standard Overlapping Panel
- Seam-to-seam, single LE protection
- Custom and Standard Sizing available

Dark gray areas indicate where panel(s) are placed to create desired protection. In this example, there is uniform 0.50mm LE protection from seam to seam (no protection changes at overlap). This is achieved through the use of additional, internal layers of lead.

Example
Order Configuration - 0.50/0.50mm LE

PROGUARD ESSENTIAL

ESSENTIAL WRAPAROUND BACK RELIEF APRON
- 0.50MM LE LEAD CONFIGURATION -



Example
Order Configuration - 0.50/0.25mm LE

ESSENTIAL

- 3/4 Overlapping Panel
- Desired protection at overlap where most sensitive organs are - helps with weight reduction
- Standard Sizing Only

Dark gray areas indicate where vest panels overlap to create desired protection. In this example, two 0.25mm vest panels overlap to create 0.50mm LE protection over the body's most sensitive organs.